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		A25	A15
GENERAL	Test name	URIC ACID	URIC ACID
	Analysis mode	endpoint mon.	endpoint mon.
	Sample type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	2	2
	Replicates	1	1
Name of assoc. constituent	-	-	
PROCEDURE	Type of reading	bichrom.	bichrom.
	Volumes		
	Sample	7.5	7.5
	Reagent 1	300	300
	Reagent 2	-	-
	Washing	1.2	1.2
	Predilution factor	-	-
	Filters		
	Main	505	505
	Reference	670	670
Times	Reading 1	300 s	312 s
	Reading 2	-	-
	Reagent 2	-	-
	Postdilution factor	2	2
CALIBRATION	Type of calibration	multiple	multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
OPTIONS	Blank absorbance limit	0.200	0.200
	Kinetic blank limit	-	-
	Linearity limit	25	25

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### BIBLIOGRAFIA

- Barham D, Trinder P. An improved colour reagent for the determination of blood glucose by oxidase system. *Analyst* 1972; 27:142-145.
- Fossati P, Prencipe L, Berti G. Use of 3,5-dichloro-2-hydroxybenzenesulfonic acid/4-aminophenazone chromogenic system in direct enzymic assay of uric acid in serum and urine. *Clin Chem* 1980; 26:227-231.
- Tietz NW. *Clinical guide to laboratory tests*, 2nd ed. Saunders Co, 1991.
- Young DS. *Effects of drugs on clinical laboratory tests*, 4th ed. AACC Press, 1995.
- Friedman and Young. *Effects of disease on clinical laboratory tests*, 3th ed. AACC Press, 1997.