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		A25	A15
GENERAL	Test name	CHOLESTEROL HDL DIRECT	CHOLESTEROL HDL DIRECT
	Analysis mode	differential bir.	differential bir.
	Sample type	serum	serum
	Units	mg/dL	mg/dL
	Reaction type	increasing	increasing
	Decimals	1	1
	Replicates	1	1
Name of assoc. constituent	-	-	
PROCEDURE	Type of reading	monoch.	monoch.
	Volumes		
	Sample	3	3
	Reagent 1	300	300
	Reagent 2	100	100
	Washing	1.2	1.2
	Predilution factor	-	-
	Filters		
	Main	535	535
	Reference	-	-
Times	Reading 1	285 s	168 s
	Reading 2	600 s	480 s
	Reagent 2	300 s	192 s
	Postdilution factor	2	2
	CALIBRATION	Type of calibration	multiple
Calibrator replicates		3	3
Blank replicates		3	3
Calibration curve		-	-

OPTIONS	Blank absorbance limit Kinetic blank limit Linearity limit	0.050 - 200	0.050 - 200
---------	--	-------------------	-------------------

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BIBLIOGRAFIA

- Warnick GR Nauck M, Rifai N. Evolution of methods for measurement of HDL-cholesterol: from ultracentrifugation to homogeneous assays. *Clin Chem* 2001; 47: 1579-96.
- National Cholesterol Education Program Expert Panel. Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (ATP III). NIH Publication. Bethesda: National Heart, Lung, and Blood Institute; 2001.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.
- Tietz NW. Clinical guide to laboratory tests, 2nd ed. Saunders Co, 1991.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.