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		A25	A15	
GENERAL	Test name	CHOLESTEROL LDL DIRECT	CHOLESTEROL LDL DIRECT	
	Analysis mode	differential bir.	differential bir.	
	Sample type	serum	serum	
	Units	mg/dL	mg/dL	
	Reaction type	increasing	increasing	
	Decimals	0	0	
	Replicates	1	1	
	Name of assoc. constituent	-	-	
	PROCEDURE	Type of reading	monoch.	monoch.
		Volumes		
Sample		3	3	
Reagent 1		300	300	
Reagent 2		100	100	
Washing		1.2	1.2	
Predilution factor		-	-	
Filters				
Main		535	535	
Reference		-	-	
Times	Reading 1	285 s	288 s	
	Reading 2	600 s	480 s	
	Reagent 2	300 s	312 s	
	Postdilution factor	2	2	

CALIBRATION	Type of calibration	multiple	multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
OPTIONS	Blank absorbance limit	0.100	0.100
	Kinetic blank limit	-	-
	Linearity limit	990	990

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BIBLIOGRAFIA

1. Nauck M, Warnick GR, Rifai N. Methods for measurement of LDL-cholesterol: a critical assessment of direct measurement by homogeneous assays versus calculation. *Clin Chem* 2002; 48: 236-54.
2. National Cholesterol Education Program Expert Panel. Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (ATP III). NIH Publication. Bethesda: National Heart, Lung, and Blood Institute; 2001.
3. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.
4. Tietz NW. Clinical guide to laboratory tests, 2nd ed. Saunders Co, 1991.
5. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.