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		A25	A15	
GENERAL	Test name	IgG	IgG	
	Analysis mode	endpoint mon.	endpoint mon.	
	Sample type	serum	serum	
	Units	mg/dL	mg/dL	
	Turbidim. test	yes	yes	
	Reaction type	increasing	increasing	
	Decimals	0	0	
Replicates	1	1		
Name of assoc. constituent	-	-		
PROCEDURE	Type of reading	monoch.	monoch.	
	Sample	3	3	
	Reagent 1	440	440	
	Reagent 2	-	-	
	Washing	1.2	1.2	
	Predilution factor	-	-	
	Filters	Main	535	535
		Reference	-	-
	Times	Reading 1	300 s	312 s
		Reading 2	-	-
Reagent 2		-	-	
Postdilution factor		1.4	1.4	
CALIBRATION	Type of calibration	multiple	multiple	
	Calibrator number	5	5	
	Calibrator replicates	3	3	
	Blank replicates	3	3	
	Calibration curve	increasing polygonal	increasing polygonal	
OPTIONS	Blank absorbance limit	0.300	0.300	
	Kinetic blank limit	-	-	
	Linearity limit	-	-	

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## BIBLIOGRAFIA

- Narayanan S. Method-comparison studies on immunoglobulins. *Clin Chem* 1982; 28: 1528-1531.
- Price CP, Spencer K and Whicher J. Light-scattering immunoassay of specific proteins: a review. *Ann Clin Biochem* 1983; 20: 1-14.
- Dati F et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference range for 14 proteins in serum based on the standardization against the IFCC/CAP reference material (CRM 470). *Eur J Clin Chem Clin Biochem* 1996; 34: 517-520.
- Young DS. Effects of drugs on clinical laboratory tests, 3th ed. AACC Press, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.