**CALIBRADOR BIOQUÍMICA**

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BIOCHEMISTRY CALIBRATOR

COD 18011 5 x 5 mL
STORE AT 2-8°C
Only for <i>in vitro</i> use in the clinical laboratory

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PRODUCT DESCRIPTION

The Biochemistry Calibrator is a lyophilized bovine serum containing component concentrations suitable for the calibration of measurement procedures.

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INSTRUCTIONS FOR USE

1. Open a vial very carefully, avoiding any loss of the lyophilized material.
2. Pipette exactly 5.00 mL of distilled water into the vial. The constituent values depend on the accuracy of this reconstitution step.
3. Close vial with the stopper and let stand for 20 minutes at room temperature.
4. Swirl gently, avoiding the formation of foam, to ensure complete dissolution of contents.
5. Use the reconstituted material as a calibrator.

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STORAGE AND STABILITY

The lyophilized Calibrator serum is stable until the expiration date given in the label.

The components of the reconstituted Calibrator serum are stable for at least 7 days at 2-8°C, excepting:

- Bilirubin is stable 8 hours at 2-8°C protected from light.
- CK is stable 6 hours at 2-8°C.
- Alkaline phosphatase values rise about 10-20% during the first 2 hours after reconstitution.
- LDH is stable 5 days at 2-8°C.

Reconstituted Calibrator Serum is stable for a least 30 days at -20°C (except for Acid Phosphatase, Bilirubin, CK and LDH that decrease their activity in at least 1-2% per week). Only freeze once.

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ASSIGNED VALUES

The assigned concentration values for components and their traceability are shown in the enclosed tables. Traceability of the results in human serum samples to reference materials or reference systems of a higher metrological order can be assured only if the BioSystems reagents and recommended measurement procedures are used.

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