## ORINA CONTROL

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## CONTROL URINE

COD 18036 1 x 20 mL	COD 18037 10 x 20 mL
Only for <i>in vitro</i> use in the clinical laboratory	

### PRODUCT DESCRIPTION

The Control Urine is a lyophilized human urine containing component concentrations suitable for the quality control of the clinical laboratories, and without preservatives which might interfere with the tests. The product is intended for intralaboratory quality control purposes only and is supplied with intervals of suggested acceptable values. The Control Urine is not to be used for demonstrating trueness of measurement.

### CONTENTS AND COMPOSITION

**Control Urine.** Lyophilized human urine (for 20 mL).

### INSTRUCTIONS FOR USE

1. Open a vial very carefully, avoiding any loss of the lyophilized material.
2. Pipette exactly 20.0 mL of distilled water into the vial. The component values depend on the accuracy of this reconstitution step.
3. Close the vial with the stopper and let stand for 20 minutes at room temperature.
4. Swirl gently, avoiding the formation of foam, to ensure complete dissolution of contents.
5. The reconstituted Control Urine is to be treated like the patient samples.

### STORAGE AND STABILITY

Store at 2-8°C.

The lyophilized Control Urine is stable until the expiration date given in the label.

The components of the reconstituted Control Urine are stable for at least 5 days at 2-8°C.

### SPECIFIED VALUES

The mean values have been obtained by interlaboratory consensus using the specified method.

The intervals of suggested acceptable values have been calculated from previous experience in interlaboratory variability and are given for orientation only; each laboratory should establish its own precision parameters.

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