SUERO CONTROL DE LÍPIDOS I

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LIPID CONTROL SERUM I

COD 18040 3 x 1 mL
STORE AT 2-8°C
Only for <i>in vitro</i> use in the clinical laboratory

PRODUCT DESCRIPTION

The Lipid Control Serum is a lyophilized defibrinated human plasma containing component concentrations suitable for the quality control of the clinical laboratories, and without preservatives which might interfere with the tests. The product is intended for intralaboratory quality control purposes only and is supplied with intervals of suggested acceptable values. The Lipid Control Serum is not to be used for demonstrating trueness of measurement.

CONTENTS AND COMPOSITION

Lipid Control Serum level I. Lyophilized defibrinated human plasma (for 1 mL).

Components from human origin have been tested and found to be negative for the presence of antibodies anti-HIV and anti-HCV, as well as for HBs antigen. However, they should be handled cautiously as potentially infectious.

INSTRUCTIONS FOR USE

1. Open a vial very carefully, avoiding any loss of the lyophilized material.
2. Pipette exactly 1.00 mL of distilled water into the vial. The component values depend on the accuracy of this reconstitution step.
3. Close the vial with the stopper and let stand for 20 minutes at room temperature.
4. Swirl gently, avoiding the formation of foam, to ensure complete dissolution of contents.
5. The reconstituted Lipid Control Serum is to be treated like the patient serum.

STORAGE AND STABILITY

The lyophilized Lipid Control Serum is stable until the expiration date given in the label.

The components of the reconstituted Lipid Control Serum are stable for at least 7 days at 2-8°C and 4 weeks at -20 °C (only freeze once).

SPECIFIED VALUES

The mean values have been obtained by interlaboratory consensus using the specified method.

The intervals of suggested acceptable values have been calculated from previous experience in interlaboratory variability and are given for orientation only; each laboratory should establish its own precision parameters.

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