PROTEIN CALIBRATORS

PROTEIN CALIBRATORS

| |
|---|
| COD 31075 5 x 1 mL |
| STORE AT 2-8°C |
| Only for <i>in vitro</i> use in the clinical laboratory |

INSTRUCTIONS FOR USE

The Calibrators are ready for use. Before opening, allow the vial to reach room temperature and then mix carefully avoiding the formation of foam.

PRODUCT DESCRIPTION

The Protein Calibrators kit consists of five liquid human serums containing protein concentrations suitable for obtaining the calibration curve necessary for the measurement of the constituents listed below.

CONTENTS AND COMPOSITION

- Calibrator level 1 with 1 mL
- Calibrator level 2 with 1 mL
- Calibrator level 3 with 1 mL
- Calibrator level 4 with 1 mL
- Calibrator level 5 with 1 mL

Components from human origin have been tested and found to be negative for the presence of antibodies anti-HIV and anti-HCV, as well as for Hbs antigen. However, they should be handled cautiously as potentially infectious.

STORAGE AND STABILITY

Store at 2-8°C.

The Protein Calibrators are stable until the expiration date printed on the label if stored tightly closed and contamination during use is avoided.

SPECIFIED VALUES

The specified mean values have been obtained by consensus from different laboratories using the indicated method. Concentration values are traceable to the Standard Reference Material 470 (Institute for Reference Materials and Measurements, IRMM).

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