**RHEUMATOID CONTROL SERUM I**

COD 31213 3 x 1 mL
STORE AT 2-8°C
Only for <i>in vitro</i> use in the clinical laboratory

PRODUCT DESCRIPTION

The Rheumatoid Control Serum is a lyophilized human serum containing rheumatoid proteins in concentrations suitable for the quality control of the clinical laboratories. No preservatives which might interfere with the tests were added. The product is intended for control of precision and accuracy of the analytical methods.

CONTENTS AND COMPOSITION

Rheumatoid Control Serum level I. Lyophilized human serum (for 1.0 mL).

Components from human origin have been tested and found to be negative for the presence of antibodies anti-HIV and anti-HCV, as well as for HBs antigen. However, they should be handled cautiously as potentially infectious.

STORAGE AND STABILITY

The lyophilized Rheumatoid Control Serum is stable until the expiration date printed on the label.

The components of the reconstituted Rheumatoid Control Serum are stable for at least 10 days at 2-8°C and 6 months at -20°C.

INSTRUCTIONS FOR USE

1. Open the vial very carefully, avoiding any loss of the lyophilized material.
2. Pipette exactly 1.0 mL of distilled water into the vial. The constituent values depend on the accuracy of this reconstitution step.
3. Replace the rubber stopper on the vial and allow the reconstituted material to stand for 20 minutes at room temperature.
4. Swirl gently. Avoid the formation of foam to ensure complete dissolution of contents.
5. Treat the Control in the analytical procedure as patient serum.

SPECIFIED VALUES

The mean value and acceptable range for each constituent have been determined by replicate analyses using the specified method. However, results may slightly differ from these values because of individual laboratory differences in technique or equipment. The acceptable ranges have been calculated from previous experience in interlaboratory variability and are given for orientation only; each laboratory should establish its own precision parameters.

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