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- Henry JB. Clinical diagnosis and management by laboratory methods. Philadelphia: WB Saunders, 1979.
 - Margolis J. The effect of colloidal silica on blood coagulation. Aust. J. Exp. Biol. 1961;39:249-58.
 - Sirridge MS. Laboratory evaluation of hemostasis. Philadelphia: Lea & Febiger, 1967.
 - Shapiro GA, Huntzinger SW, Wilson JE. Variation among commercial activated partial thromboplastin time reagents in response to heparin. Am. J. Clin. Pathol. 1977;67:477.
 - Van den Besselaar AMHP, Neuteboom J, and Bertina RM. Effect of Synthetic Phospholipids on the response of the Activated Partial Thromboplastin Time to Heparin. Blood Coag. Fibrin. 1993; 4: 895-903.
 - Ray MJ, and Hawson GAT. A Comparison of Two APTT Reagents Which Use Silica Activators. Clin. Lab. Haemat. 1989; 11:221-232.
 - Turi DC, and Peerschke EI. Sensitivity of Three Activated Partial Thromboplastin Time Reagents to Coagulation Factor Deficiencies. Am. J. Clin. Pathol. 1986; 85:43-49.
 - Van der Velde EA, and Poller L. The APTT Monitoring of Heparin - The ISTH/ICSH Collaborative Study. Thromb. Haemos. 1995; 73:73-81.
 - Brandt JT, Triplett DA, Rock WA, Bovill EG, and Arkin CF. Effect of Lupus Anticoagulants on the Activated Partial Thromboplastin Time. Arch. Pathol. Lab. Med. 1991; 115:109-114.
 - Clinical and Laboratory Standards Institute/NCCLS. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays, Fourth Edition, NCCLS/CLSI Document H21-A4; Vol. 23 No. 35.
 - Clinical and Laboratory Standards Institute/NCCLS. One Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test; Approved Guideline, NCCLS/CLSI Document H47-A; Vol. 16 No. 3.
 - Zucker S, Cathey MH, and West B. Preparation of Quality Control Specimen for Coagulation. Am. J. Clin. Pathol. 1970; 53: 924-927.
 - Westgard JO, and Barry PL. Cost-Effective Quality Control; Managing the Quality and Productivity of Analytical Process, AACC Press 1986.
 - Hirsh J, Raschke R, Warkentin TE, Dalen JE, Deykin D, Poller L. Heparin: mechanism of action, pharmacokinetics, dosing considerations, monitoring, efficacy, and safety. Chest 1995 Oct; 108 (4 Suppl): 258S-275S.
 - Brill-Edwards P, Ginsburg JS, Johnson M, and Hirsch J. Establishing a Therapeutic Range for Heparin Therapy. Ann. Int. Med. 1993; 119: 104-109.
- ACL, ACL Futura, ACL TOP and ELECTRA are registered trademarks of Instrumentation Laboratory.
- 2002 Instrumentation Laboratory
Issued June 2006

Symbols used / Verwendete Symbole / Símbolos utilizados / Symboles utilisés / Simboli impiegati / Símbolos utilizados / Anvendte symboler / Använda Symboler / Χρησιμοποιηθέντα σύμβολα

<p>IVD</p> <p><i>In vitro</i> diagnostic medical device <i>In-vitro</i> Diagnostikum De uso diagnóstico <i>in vitro</i> Dispositif médical de diagnostic <i>in vitro</i> Per uso diagnostico <i>in vitro</i> Dispositivo médico para utilização em diagnóstico <i>in vitro</i> "in vitro" diagnostisk udstyr <i>In vitro</i> diagnostisk medicinsk produkt Προϊόν για διαγνωστική χρήση <i>In vitro</i></p>	<p>LOT</p> <p>Batch code Chargen-Bezeichnung Identificación número de lote Désignation du lot Numero del lotto Número de lote Batch nr. Tillverkningskod Αρ. Παρτίδας</p>	<p></p> <p>Use by Verwendbar bis Caducidad Utilisable jusqu'à Da utilizzare prima del Data límite de utilização Anvendelse Användning Χρήση έως</p>	<p></p> <p>Temperature limitation Festgelegte Temperatur Temperatura de Almacenamiento Températures limites de conservation Limiti di temperatura Limite de temperatura Temperatur begrænsninger Temperatur gräns Περιορισμοί θερμοκρασίας</p>	<p></p> <p>Consult instructions for use Beilage beachten Consultar la metódica Lire le mode d'emploi Vedere istruzioni per l'uso Consultar as instruções de utilização Se vejledning for anvendelse Ta del av instruktionerna före användning Συμβουλευτείτε τις οδηγίες χρήσης</p>	<p>CONTROL</p> <p>Control Kontrollen Control Contrôle Controllo Controlo Kontrol Kontroll Υλικό ποιοτικού ελέγχου</p>	<p></p> <p>Biological risks Biologisches Risiko Riesgo biológico Risque biologique Rischio biologico Risco biológico Miljø oplysninger Biologiska risker Βιολογικοί κίνδυνοι</p>	<p></p> <p>Manufacturer Hergestellt von Fabricado por Fabricant Prodotto da Fabricado por Producent Tillverkare Κατασκευαστής</p>	<p>EC REP</p> <p>Authorised representative Bevollmächtigter Representante autorizado Mandataire Rappresentanza autorizzata Representante autorizado Leverandør Auktoriserad representant Εξουσιοδοτημένος αντιπρόσωπος</p>
--	--	--	---	---	--	---	--	---