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<p>1. Van den Besselaar AMHP, Neuteboom J, and Bertina RM. Effect of Synthetic Phospholipids on the response of the Activated Partial Thromboplastin Time to Heparin. <i>Blood Coag. Fibrin.</i> 1993; 4: 895-903.</p> <p>2. Ray MJ, and Hawson GAT. A Comparison of Two APTT Reagents Which Use Silica Activators. <i>Clin. Lab. Haemat.</i> 1989; 11:221-232.</p> <p>3. Turi DC, and Peerschke EI. Sensitivity of Three Activated Partial Thromboplastin Time Reagents to Coagulation Factor Deficiencies. <i>Am J. Clin. Pathol.</i> 1986; 85:43-49.</p> <p>4. Van der Velde EA, and Poller L. The APTT Monitoring of Heparin - The ISTH/ICSH Collaborative Study. <i>Thromb. Haemos.</i> 1995; 73:73-81.</p> <p>5. Brandt JT, Triplett DA, Rock WA, Bovill EG, and Arkin CF. Effect of Lupus Anticoagulants on the Activated Partial Thromboplastin Time. <i>Arch. Pathol. Lab. Med.</i> 1991; 115:109-114.</p>	<p>6. National Committee for Clinical Laboratory Standards. <i>Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays</i>, Fourth Edition, NCCLS Document H21-A4; Vol. 23 No. 35.</p> <p>7. National Committee for Clinical Laboratory Standards. <i>One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test; Approved Guideline</i>, NCCLS Documents H47-A, Vol. 16 No. 3.</p> <p>8. Zucker S, Cathey MH, and West B. Preparation of Quality Control Specimen for Coagulation. <i>Am. J. Clin. Pathol.</i> 1970; 53:924-927.</p> <p>9. Westgard JO, and Barry PL. <i>Cost-Effective Quality Control: Managing the Quality and Productivity of Analytical Process</i>. AACC Press 1986.</p>	<p>10. Solberg HE. Approved Recommendation (1987) on the Theory of Reference Values. Part 5. Statistical Treatment of Collected Reference Values. Determination of Reference Limits. <i>J Clin Chem Clin Biochem</i> 1987; 25: 645-656.</p> <p>11. Hirsh J, Raschke R, Warkentin TE, Dalen JE, Deykin D, Poller L. Heparin: mechanism of action, pharmacokinetics, dosing considerations, monitoring, efficacy, and safety. <i>Chest</i> 1995 Oct; 108 (4 Suppl): 258S-275S.</p> <p>12. Brill-Edwards P, Ginsburg JS, Johnson M, and Hirsch J. Establishing a Therapeutic Range for Heparin Therapy. <i>Ann. Int. Med.</i> 1993; 119:104-109.</p>
---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

ACL, ACL Futura and ACL TOP are trademarks of Instrumentation Laboratory

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Printed in U.S.A.