Thrombin Time - 0009758515

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1. Biggs R, Rizza C. Human Blood Coagulation, Haemostasis and Thrombosis. III ed. Oxford, England: Blackwell Scientific Publications, 1984.
2. Colman R, *et al.* Hemostasis and Thrombosis, Basic Principles and Clinical Practice. Philadelphia: J.B. Lippincott Co., 1987.
3. National Committee for Clinical Laboratory Standards. Preparation and Testing of Reagent Water in the Clinical Laboratory, Third Edition, NCCLS Document C3-A3; Vol. 17 No. 18.
4. National Committee for Clinical Laboratory Standards. Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays, Fourth Edition, NCCLS Document H21-A4; Vol. 23 No. 35.
5. Zucker S, Cathey MH, West B. Preparation of Quality Control Specimens for Coagulation, Am. J. Clin. Pathol. 1970; 53: 924-927.
6. Westgard JO, and Barry PL. Cost-Effective Quality Control: Managing the Quality and Productivity of Analytical Process, AACC Press 1986.
7. IFCC Committee Report, J. Clin. Chem. Clin. Biochem. 1983; 21: 749-760.

ACL, ACL Futura and ACL TOP are trademarks of Instrumentation Laboratory.

©1998 Instrumentation Laboratory
Issued March 2005
Printed in U.S.A.

Symbols used / Verwendete Symbole / Símbolos utilizados / Symboles utilisés / Simboli impiegati / Símbolos utilizados / Anvendte symboler / Använda Symboler / Χρησιμοποιηθέντα σύμβολα

<p>IVD</p> <p><i>In vitro</i> diagnostic medical device In-vitro Diagnostikum De uso diagnóstico <i>in vitro</i> Dispositif médical de diagnostic <i>in vitro</i> Per uso diagnostico <i>in vitro</i> Dispositivo médico para utilização em diagnóstico <i>in vitro</i> "in vitro" diagnostisk udstyr In vitro diagnostisk medicinsk produkt</p> <p>Προϊόν για διαγνωστική χρήση <i>In vitro</i></p>	<p>LOT</p> <p>Batch code Chargen-Bezeichnung Identificación número de lote Désignation du lot Numero del lotto Número de lote Batch nr. Tillverkningskod Αρ. Παρτίδας</p>	<p></p> <p>Use by Verwendbar bis Caducidad Utilisable jusqu'à Da utilizzare prima del Data limite de utilização Anvendelse Användning Χρήση έως</p>	<p></p> <p>Temperature limitation Festgelegte Temperatur Temperatura de Almacenamiento Températures limites de conservation Limiti di temperatura Limite de temperatura Temperatur begrænsninger Temperatur gräns Περιορισμοί θερμοκρασίας</p>	<p></p> <p>Consult instructions for use Beilage beachten Consultar la metódica Lire le mode d'emploi Vedere istruzioni per l'uso Consultar as instruções de utilização Se vejledning for anvendelse Ta del av instruktionen före användning Συμβουλευτείτε τις οδηγίες χρήσης</p>	<p>CONTROL</p> <p>Control Kontrollen Control Contrôle Controllo Controllo Kontrol Kontroll Υλικό ποιοτικού ελέγχου</p>	<p></p> <p>Biological risks Biologisches Risiko Riesgo biológico Risque biologique Rischio biologico Risco biológico Miljø oplysninger Biologiska risker Βιολογικοί κίνδυνοι</p>	<p></p> <p>Manufacturer Hergestellt von Fabricado por Fabricant Prodotto da Fabricado por Producent Tillverkare Κατασκευαστής</p>	<p>EC REP</p> <p>Authorised representative Bevollmächtigter Representante autorizado Mandataire Rappresentanza autorizzata Representante autorizado Leverandør Auktoriserad representant Εξουσιοδοτημένος αντιπρόσωπος</p>
---	--	--	---	--	---	---	--	---