# Antithrombin - 0020008900

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- Wells PS, Blajchman MA, Henderson P. The prevalence of Antithrombin Deficiency in Healthy Blood Donors. A Cross-Sectional Study, *Am. J. Hematol.* 1994; 4: 321-324.
- Tait RC, Walker ID, Perry DJ, Carrell RW, Islam SIA, McCall F, Mitchell R, Davidson JF. Prevalence of Antithrombin III Deficiency Subtypes in 4000 Healthy Blood Donors, *Thromb. Haemostas.* 1991; 65: 534 (abst).
- Harper PL, Luddington RJ, Daly M. The Incidence of Dysfunctional Antithrombin Variants: Four Cases in 210 Patients with Thromboembolic Disease, *Br. J. Haematol.* 1991; 77: 360-364.
- Perry DJ. Antithrombin and its Inherited Deficiencies, *Blood.* 1994; 8: 37-55.
- Pabinger I, Brucker S, Kyrle PA. Hereditary Deficiency of Antithrombin III, Protein C and Protein S. Prevalence in Patients with a History of Venous Thrombosis and Criteria for Rational Patient Screening, *Blood Coag. Fibrinolys.* 1992; 3: 547-553.
- Fourrier F, Chopin C, Goudemand J, Hendrycx S, Caron C, Rime A, Marey A, Lestavel P. Septic Shock, Multiorgan Failure, and Disseminated Intravascular Coagulation. Compared Patterns of Antithrombin III, Protein C and Protein S Deficiencies, *Chest.* 1993; 101: 816-823.
- Knot E, Ten Cate JW, Drijfhout HR, Kahlé LH, Tytgat GN. Antithrombin III Metabolism in Patients with Liver Disease, *J. Clin. Pathol.* 1984; 37: 523-530.
- Schulman S, Tengbom L. Treatment of Venous Thromboembolism, in Patients with Congenital Deficiency of Antithrombin, *Thromb. Haemostas.* 1992; 68: 634-636.
- Schwartz RS, Bauer KA, Rosenberg RD, Kavanaugh EJ, Davies DC, Bogdanoff DA. Clinical Experience with Antithrombin III Concentrate in treatment of Congenital and Acquired Deficiency of Antithrombin, *Am. J. Med.* 1989; 87 (sup 3B): 53-60.
- Bohner J, von Pape K, Blaurock M. Thrombin Based Antithrombin Assays show Over Estimation of Antithrombin III Activity in Patients on Heparin Therapy due to Heparin Cofactor II Influence, *Thromb. Haemostas.* 1994; 71: 280-283.
- Demers C, Henderson P, Blajchman MA, Wells MJ, Mitchell L, Johnston M, Ofusu FA, Fernandez-Rachubinski F, Andrew M, Hirsh J, Ginsberg J. An Antithrombin III Assay Based on Factor Xa Inhibition Provides a more Reliable Test to Identify Congenital Antithrombin III Deficiency than an Assay Based on Thrombin Inhibition, *Thromb. Haemostas.* 1993; 69: 231-235.
- Clinical and Laboratory Standards Institute/NCCLS. Preparation and Testing of Reagent Water in the Clinical Laboratory, Third Edition, CLSI/NCCLS Document C3-A3; Vol. 17 No. 18.
- Clinical and Laboratory Standards Institute/NCCLS. Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays, Fourth Edition, CLSI/NCCLS Document H21-A4; Vol. 23 No. 35.
- Zucker S, Cathey MH, West B. Preparation of Quality Control Specimens for Coagulation, *Am J. Clin. Pathol.* 1970; 53: 924-927.
- Westgard JO and Barry PL. Cost-Effective Quality Control: managing the Quality and Productivity of Analytical Process, AACC Press 1986.
- Kottke-Marchant K, Duncan A. Antithrombin Deficiency: Issues in Laboratory Diagnosis, *Arch Pathol Lab Med.* 2002; 126:1326-1336.

ACL and ACL FUTURA are trademarks of Instrumentation Laboratory

©1999 Instrumentation Laboratory  
Issued January 2007

## Symbols used / Verwendete Symbole / Símbolos utilizados / Symboles utilisés / Simboli impiegati / Símbolos utilizados

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