Liquid Antithrombin - 0020002500

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- Wells PS, Blajchman MA, Henderson P. The prevalence of Antithrombin Deficiency in Healthy Blood Donors. A Cross-Sectional Study, Am. J. Hematol. 1994; 4: 321-324.
- Tait RC, Walker ID, Perry DJ, Carrell RW, Islam SIA, McCall F, Mitchell R, Davidson JF. Prevalence of Antithrombin III Deficiency Subtypes in 4000 Healthy Blood Donors, Thromb. Haemost. 1991; 65: 534 (abst).
- Harper PL, Luddington RJ, Daly M. The Incidence of Dysfunctional Antithrombin Variants: Four Cases in 210 Patients with Thromboembolic Disease, Br. J. Haematol. 1991; 77: 360-364.
- Perry DJ. Antithrombin and its Inherited Deficiencies, Blood. 1994; 8: 37-55.
- Pabinger I, Brucker S, Kyrle PA. Hereditary Deficiency of Antithrombin III, Protein C and Protein S. Prevalence in Patients with a History of Venous Thrombosis and Criteria for Rational Patient Screening, Blood Coag. Fibrinolys. 1992; 3: 547-553.
- Fourrier F, Chopin C, Goudemand J, Hendrycx S, Caron C, Rime A, Marey A, Lestavel P. Septic Shock, Multiorgan Failure, and Disseminated Intravascular Coagulation. Compared Patterns of Antithrombin III, Protein C and Protein S Deficiencies, Chest. 1993; 101: 816-823.
- Knot E, Ten Cate JW, Drijfhout HR, Kahlé LH, Tytgat GN. Antithrombin III Metabolism in Patients with Liver Disease, J. Clin. Pathol. 1984; 37: 523-530.
- Schulman S, Tengborn L. Treatment of Venous Thromboembolism, in Patients with Congenital Deficiency of Antithrombin, Thromb. Haemostas. 1992; 68: 634-636.
- Schwartz RS, Bauer KA, Rosenberg RD, Kavanaugh EJ, Davies DC, Bogdanoff DA. Clinical Experience with Antithrombin III Concentrate in treatment of Congenital and Acquired Deficiency of Antithrombin, Am. J. Med. 1989; 87 (sup 3B): 53-60.
- Bohner J, von Pape K, Blaurock M. Thrombin Based Antithrombin Assays show Over Estimation of Antithrombin III Activity in Patients on Heparin Therapy due to Heparin Cofactor II Influence, Thromb. Haemost. 1994; 71: 280-283.
- Demers C, Henderson P, Blajchman MA, Wells MJ, Mitchell L, Johnston M, Ofusu FA, Fernandez-Rachubinski F, Andrew M, Hirsh J, Ginsberg J. An Antithrombin III Assay Based on Factor Xa Inhibition Provides a more Reliable Test to Identify Congenital Antithrombin III Deficiency than an Assay Based on Thrombin Inhibition, Thromb. Haemost. 1993; 69: 231-235.
- Clinical and Laboratory Standards Institute/NCCLS. Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays, Fourth Edition, CLSI/NCCLS Document H21-A4; Vol. 23 No. 35.
- Zucker S, Cathey MH, West B. Preparation of Quality Control Specimens for Coagulation, Am J. Clin. Pathol. 1970; 53: 924-927.
- Westgard JO and Barry PL. Cost-Effective Quality Control: managing the Quality and Productivity of Analytical Process, AACC Press 1986.
- Kottke-Marchant K, Duncan A. Antithrombin Deficiency: Issues in Laboratory Diagnosis, Arch Pathol Lab Med. 2002; 126:1326-1336.

ACL, ACL Futura and ACL TOP™ are trademarks of Instrumentation Laboratory

©2000 Instrumentation Laboratory
Issued October 2006

Symbols used / Verwendete Symbole / Símbolos utilizados / Symboles utilisés / Simboli impiegati / Símbolos utilizados / Anvendte symboler / Använda Symboler / Χρησιμοποιηθέντα σύμβολα

IVD	LOT				CONTROL			EC REP
<i>In vitro</i> diagnostic medical device <i>In vitro</i> Diagnostikum De uso diagnóstico <i>in vitro</i> Dispositif médical de diagnostic <i>in vitro</i> Per uso diagnostico <i>in vitro</i> Dispositivo médico para utilização em diagnóstico <i>in vitro</i> "in vitro" diagnostisk udstyr <i>In vitro</i> diagnostisk medicinsk produkt Προϊόν για διαγνωστική χρήση <i>In vitro</i>	Batch code Chargen-Bezeichnung Identificación número de lote Désignation du lot Numero del lotto Número de lote Batch nr. Tilverkningskod Αρ. Παρτίδας	Use by Verwendbar bis Caducidad Utilisable jusqu'à Da utilizzare prima del Data limite de utilização Anvendelse Användning Χρήση έως	Temperature limitation Festgelegte Temperatur Temperatura de Almacenamiento Températures limites de conservation Limiti di temperatura Límite de temperatura Temperatur begrænsninger Temperatur gräns Περιορισμοί θερμοκρασίας	Consult instructions for use Beilage beachten Consultar la metódica Lire le mode d'emploi Vedere istruzioni per l'uso Consultar as instruções de utilização Se vejledning for anvendelse Ta del av instruktionerna före användning Συμβουλευτήτε τις οδηγίες χρήσης	Control Kontrollen Control Contrôle Controllo Controllo Kontrol Κontrol Υλικό ποιοτικού ελέγχου	Biological risks Biologisches Risiko Riesgo biológico Risque biologique Rischio biologico Risco biológico Miljø oplysninger Biologiska risker Βιολογικοί κίνδυνοι	Manufacturer Hergestellt von Fabricado por Fabricant Prodotto da Fabricado por Producent Tilverkare Κατασκευαστής	Authorised representative Bevollmächtigter Representante autorizado Mandataire Rappresentanza autorizzata Representante autorizado Leverandør Auktoriserad representant Εξουσιοδοτημένος αντιπρόσωπος