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1. Love PE, Santoro SA. Antiphospholipid Antibodies: Anti Cardiolipin and the Lupus Anticoagulant in Systemic Lupus Erythematosus (SLE) and in Non-SLE Disorders, Ann. Int. Med. 1990; 112: 682-698.
2. Exner T, Triplett DA, Taberner D, Machin SJ. Guidelines for Testing and Revised Criteria for Lupus Anticoagulants, Thromb. Haemostas. 1991; 65: 320-322.
3. Rauch J, Tannebaum M and Janoff AS. Distinguishing Plasma Lupus Antifactor Antibodies using Hexagonal Phase II Phospholipids, Thromb. Haemostas. 1989; 62: 892-896.
4. Feinstein DI. Lupus Anticoagulant, Thrombosis and Fetal Loss, N. Eng. J. Med. 1985; 313: 1248-1350.
5. Thiagarajan P, Pengo V, Shapiro SS. The use of the Dilute Russell's Viper Venom Time for the Diagnosis of Lupus Anticoagulants, Blood Coag. Fibrinol. 1990; 1: 259-266.
6. National Committee for Clinical Laboratory Standards. Preparation and Testing of Reagent Water in the Clinical Laboratory, Third Edition, NCCLS Document C3-A3; Vol. 17 No. 18.
7. National Committee for Clinical Laboratory Standards. Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays, Fourth Edition, NCCLS Document H21-A4; Vol. 23 No. 35.
8. Zucker S, Cathey MH, West B. Preparation of Quality Control Specimen for Coagulation, Am. J. Pathol. 1970; 53: 924-927.
9. Westgard JO and Barry PL. Cost-Effective Quality Control: Managing the Quality and Productivity of Analytical Process, AACC Press 1986.
10. Brandt JT, Triplett DA, Alving B, Scharrer I, Criteria for the Diagnosis of Lupus Anticoagulants: An Update. Thromb.Haemostas. 1995, 74 (4), 1185-90

ACL, ACL Futura and ACL TOP are trademarks of Instrumentation Laboratory Issued February 2005

Symbols used / Verwendete Symbole / Símbolos utilizados / Symboles utilisés / Simboli impiegati / Símbolos utilizados / Anvendte symboler / Använda Symboler / Χρησιμοποιηθέντα σύμβολα

<p><b>IVD</b></p> <p><i>In vitro</i> diagnostic medical device  <i>In-vitro</i> Diagnostikum          De uso diagnóstico <i>in vitro</i>          Dispositif médical de diagnostic <i>in vitro</i>          Per uso diagnóstico <i>in vitro</i>          Dispositivo médico para utilização <i>in vitro</i>          "in vitro" diagnostisk udstyr  <i>In vitro</i> diagnostisk medicinsk produkt</p> <p>Προϊόν <del>in vitro</del> γνωστικής χρήσης</p>	<p><b>LOT</b></p> <p>Batch code          Chargen-Bezeichnung          Caducidad          Identificación número de lote          Désignation du lot          Numero del lotto          Número de lote          Batch nr.          Tillverkningskod</p> <p>Αρ. Παρτίδας</p>	<p></p> <p>Use by          Verwendbar bis          Caducidad          Utilisable jusqu'à          Da utilizzare prima del          Data limite de utilização          Anvendelse          Användning</p> <p>Χρήση έως</p>	<p></p> <p>Temperature limitation          Festgelegte          Temperatur          Temperatura de          Almacenamiento          Températures limites de conservation          Limiti di temperatura          Limite de temperatura          Temperatur          begrænsninger          Temperatur gråns</p> <p>Περιορισμοί θερμοκρασίας</p>	<p></p> <p>Consult instructions for use          Beilage beachten          Consultar la metódica          Lire le mode d'emploi          Vedere istruzioni per l'uso          Consultar as instruções de utilização          Se vejledning for anvendelse          Ta del av instruktionen före användning</p> <p>Συμβουλευτείτε τις οδηγίες χρήσης</p>	<p><b>CONTROL</b></p> <p>Control          Kontrollen          Control          Contrôle          Controllo          Controllo          Kontrol          Kontroll</p> <p>Υλικό ποιοτικού ελέγχου</p>	<p></p> <p>Biological risks          Biologisches Risiko          Riesgo biológico          Risque biologique          Rischio biologico          Risco biológico          Miljø oplysninger          Biologiska risker</p> <p>Βιολογικοί κίνδυνοι</p>	<p></p> <p>Manufacturer          Hergestellt von          Fabricado por          Fabricant          Prodotto da          Fabricado por          Producent          Tillverkare</p> <p>Κατασκευαστής</p>	<p><b>EC REP</b></p> <p>Authorised representative          Bevollmächtigter          Representante autorizado          Mandataire          Rappresentanza autorizzata          Representante autorizado          Leverandør          Auktoriserad representant</p> <p>Εξουσιοδοτημένος αντιπρόσωπος</p>
--	---	--	--	--	---	---	--	---