Coagulation

IL Test™ von Willebrand Factor 20002300

INSERT REVISION 09/2001

2 vials x 3 mL Latex Reagent
2 vials x 4 mL Reaction Buffer

Intended use

Automated latex enhanced immunoassay for the quantitative determination of von Willebrand Factor Antigen (VWF:Ag) in human citrated plasma on IL Coagulation Systems.

Summary and principle

The diagnosis of von Willebrand disease (VWD), probably the most common congenital bleeding disorder, requires a number of special testing at the laboratory level.¹ Among them, VWF:Ag determination is essential and must be performed on every patient to reach a proper diagnosis.¹ Depending upon these laboratory findings, VWD is classified into type 1 (the most frequent form being 70-80% of VWD), type 2 or type 3 (1 to 3% of VWD) groups.^{2,3} Type 1 shows a reduction of VWF although its structure and functionality is normal. In type 3, VWF is almost absent in plasma. In type 2 the quantity of VWF in plasma may be normal or slightly reduced but its molecular structure and its functionality is abnormal. Type 2 may be further characterized into subtypes by multimeric structure analysis of VWF. Apart from the above described inherited VWD, acquired VWD due to autoantibodies or to various disease states resulting in low rates of VWF synthesis has been reported.

On the other side, chronic or acute inflammatory diseases or processes involving damage of the vascular endothelium yield abnormally high concentrations of VWF.⁴

IL Test™ VWF:Ag is a latex particle enhanced immunoturbidimetric assay to quantify VWF:Ag in plasma. When a plasma containing VWF:Ag is mixed with the Latex Reagent and the Reaction Buffer included in the kit, the coated latex particles agglutinate. The degree of agglutination is directly proportional to the concentration of VWF:Ag in the sample and is determined by measuring the decrease of transmitted light caused by the aggregates.⁵

Composition

The IL VWF:Ag kit consists of:

Latex Reagent: 2 vials of a suspension of polystyrene latex particles coated with a rabbit polyclonal antibody directed against VWF containing bovine serum albumin, buffer, stabilizer and preservative.

Reaction Buffer: 2 vials of HEPES buffer containing bovine serum albumin, stabilizers and preservative.

CAUTION: All reagents contain less than 0.1% sodium azide that may form explosive azides in metal plumbing. Use proper disposal procedures.

This product is for *in vitro* diagnostic use.

Preparation

Latex Reagent: Invert to mix before use.

Reaction Buffer: Invert to mix before use.

Reagent storage and stability

Unopened reagents are stable until the expiration date shown on the vial when stored at 2-8°C.

Opened reagents are stable 3 months at 2-8°C in the original vial or 1 week at 15°C on the ACL 9000 and on ACL Futura/ACL Advance Systems. Do not freeze.

For optimal stability remove reagents from the system and store them at 2-8°C in the original vial.

Instrument/test procedures

Refer to the appropriate IL instrument's Operator's Manual and/or Application Manual for the complete assay procedure instructions.

Specimen collection and preparation

Nine parts of freshly drawn venous blood are collected into one part trisodium citrate (3.2%). Refer to NCCLS Document H21-A3 for further instructions on specimen collection, handling and storage.⁶ Thaw frozen specimens at 37°C for at least 15 minutes and centrifuge plasma before testing. After thawing the assay must be performed within 2 hours.

Additional reagents and control plasmas

The following are not supplied with the kit and must be purchased separately.

	Americas and Pacific Rim Cat. No.	Europe Cat. No.
Calibration plasma	20000000	08467300
Abnormal chromogenic control plasma Level 1/2	08467600	08467600
Factor diluent	09757600	09757600

Quality control

Two levels of control are recommended for a complete quality control program.⁷ Abnormal chromogenic control plasma Level 1 (that contains a level of VWF:Ag within the normal range) and Abnormal chromogenic control plasma Level 2 (that contains a level of VWF:Ag within the abnormal range) are designed for this program. Each laboratory should establish its own mean and standard deviation and should establish a quality control program to monitor laboratory testing. Controls should be analyzed at least once every 8 hour shift in accordance with good laboratory practice. Refer to the instrument's Operator's Manual for additional information. Refer to Westgard *et al* for identification and resolution of out-of-control situations.⁸

Results

VWF:Ag results are reported in % of normality.

Refer to the instrument's Operator's Manual for additional information. The assay results should be used with other information, including the clinical context, in forming a diagnosis.

Limitations/interfering substances

VWF:Ag results are not affected by hemoglobin up to 80 mg/dL, bilirubin up to 10 mg/dL and lipids up to 1000 mg/dL. The presence of Rheumatoid Factor may produce an overestimation of the test result. Hemolyzed samples should not be assayed.

Expected values

A normal range study was performed using IL VWF:Ag kit on IL Coagulation Systems.

Blood ABO Type ⁹	N	ACL % VWF:Ag	ACL Futura/ACL Advance % VWF:Ag
O	58	41.1 - 125.9	40.6 - 122.0
A + B + AB	56	61.3 - 157.8	61.3 - 152.6

These results were obtained using several lots of reagents. Ranges (mean ± 2SD) were estimated testing plasma samples from healthy adult blood bank donors. Due to many variables which may affect results, each laboratory should establish its own VWF:Ag normal range.

Performance characteristics

Precision:

Within run precision was assessed over multiple runs.

ACL

Mean (% VWF:Ag)	33.7	80.8	101.5
CV%	2.2	1.3	1.4

ACL Futura/ACL Advance

Mean (% VWF:Ag)	35.5	79.5	99.5
CV%	3.2	2.5	3.5

Correlation:

System	slope	intercept	r	Comparative method
ACL	1.00	-2.2	0.996	VWF:Ag EIA
ACL Futura/ACL Advance	1.03	-3.8	0.997	VWF:Ag EIA

The precision and correlation results were obtained using specific lots of reagents and controls.

Linearity:

System

ACL and ACL Futura/ACL Advance 10 - 150 (% VWF:Ag)

Without prozone effect up to 1600%.

Sample results above 150% should be diluted 1:4 with Factor diluent (100 µL of sample + 300 µL of Factor diluent) and reassayed. This procedure should be repeated if the result is again above 600% (dilution ratio 1:16). The printed results must be multiplied by 4 or 16 (depending upon the number of dilution steps performed) to correct for the dilution. Instruments with Auto Rerun capability perform the 1:4 dilution automatically expanding the test range to 600%. Above 600%, samples must be manually diluted 1:16, reassayed and results multiplied by the corresponding factor.

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1. Montgomery RR and Collier BS. Von Willebrand Disease. In: Hemostasis and Thrombosis: Basic Principles and Clinical Practice. Third edition. Chapter 7. Edited by Colman RW, Hirsh J, Marder VJ and Salzman EW. JB Lippincott Company, Philadelphia. 1994.
2. Ruggeri ZM. Structure and Function of von Willebrand Factor: Relationship to von Willebrand's Disease, Mayo Clin. Proc. 1991; 66: 847-861.
3. Sadler JE. A Revised Classification of von Willebrand Disease, Thromb. and Haemostas. 1994; 71: 520-525.
4. Mannucci PM. von Willebrand Factor - A Marker of Endothelial Damage? Arterioscler. Thromb. Vasc. Biol. 1998; 18: 1359-1362.
5. Newman DJ, Henneberry H and Price CP. Particle enhanced light scattering immunoassay, Ann. Clin. Biochem. 1992; 29: 22-42.
6. National Committee for Clinical Laboratory Standards. Collection, Transport and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays; Approved Guidelines, Third Edition, NCCLS Document H21-A3; Vol. 18 No. 20.
7. Zucker S, Cathey MH and West B. Preparation of Quality Control Specimens for Coagulation, Am. J. Clin. Pathol. 1970; 53: 924-927.
8. Westgard JO and Barry PL. Cost-Effective Quality Control: Managing the Quality and Productivity of Analytical Processes, AACC Press 1986.
9. Gill JC, Endres-Brooks J, Bauer PJ, Marks WJ and Montgomery RR. The Effect of ABO Blood Group on the Diagnosis of von Willebrand Disease, Blood. 1987; 69: 1691-1695.

Tests on the ACL Systems

VWF:Ag 50 tests (approx.)

Tests on the ACL Futura/ACL Advance Systems

VWF:Ag 60 tests (approx.)

IL Test, ACL and ACL Futura/ACL Advance are trademarks of Instrumentation Laboratory.
©1999 Instrumentation Laboratory.
Issued September 2001.
Printed in Spain.
3800-2723 R03 09/2001