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1. Refsum H, Smith AD, Ueland PM, Nexø E, Clarke R, McPartlin J, Johnston C, Engbaek F, Schneede J, McPartlin C, Scott JM. Facts and recommendations about total homocysteine determinations: an expert opinion. Clin Chem 2004; 50 (1):3-32
2. Finkelstein JD. Methionine metabolism in mammals. J Nutr Biochem 1990; 1:228-237.
3. Mudd SH, et al. Disorder of Transsulfuration. In: C.R. Scriver, A.L. Beaudet, W.S. Sly, et al., eds. The metabolic basis of inherited disease New York: McGraw-Hill, 1995: 1297-1327.
4. Boushey, Carol J, et al. A Quantitative Assessment of Plasma Homocysteine as a Risk Factor for Vascular Disease, JAMA, October 4, 1995-Vol 274, No. 13: 1049-1057.
5. Allen, Robert H, et al. Diagnosis of Cobalamin Deficiency I: Usefulness of Serum Methylmalonic Acid and Total Homocysteine Concentrations, American Journal of Hematology 1990; 34:90-98.
6. Stabler, Sally P, et al. The Use of Homocysteine and Other Metabolites in the Specific Diagnosis of Vitamin B-12 Deficiency J Nutr 126: 1266S-1272S, 1996.
7. Guttormsen AB, Svarstad E, Ueland PM, Refsum H. Elimination of homocysteine from plasma in subjects with end-stage renal failure. Irish J Med Sci 1995; 164 (Suppl. 15): 8.
8. Bostom AG, Lathrop L. Hyperhomocysteinemia in end-stage renal disease (ESDR): Prevalence, etiology, and potential relationship to arteriosclerotic outcomes. Kidney Int 1997; 52:10-20.
9. Clinical Laboratory Standards Institute/NCCLS. Preparation and Testing of Reagent Water in the Clinical Laboratory, Third Edition, CLSI/NCCLS Document C3-A3; Vol. 17 No. 18.
10. Clinical Laboratory Standards Institute/NCCLS. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays, Fourth Edition, CLSI/NCCLS Document H21-A4; Vol. 23 No. 35.
11. Zucker S, Cathey MH, West B. Preparation of Quality Control Specimens for Coagulation, Am J Clin Pathol 1970; 53: 924-927.
12. Westgard JO, Barry PL. Cost-Effective Quality Control: Managing the Quality and Productivity of Analytical Processes, AACC Press 1986.
13. Refsum H, Ueland PM, Svardal M. Fully Automated Fluorescence Assay for Determining Total Homocysteine in Plasma. Clin Chem 1989; 35:1921-27.
14. Fiskerstrand T, Refsum H, Kvalheim G, Ueland PM. Homocysteine and other thiols in plasma and urine: Automated determination and sample stability. Clin Chem 1993; 39: 263-271.
15. Primus FJ, Kelly EA, Hansen HJ, et al. "Sandwich"- Type Immunoassay of Carcinoembryonic Antigen in Patients Receiving Murine Monoclonal Antibodies for Diagnosis and Therapy. J Clin Chem 1988; 34: 261-4.
16. Schroff RW, Foon KA, Beatty SM et al. Human Anti-Murine Immunoglobulin Responses in Patients receiving Monoclonal Antibody Therapy. Cancer Res 1985; 45: 879-85.
17. Davis RE. Clinical chemistry of vitamin B12. Adv Clin Chem 1985; 24: 163-216.
18. Compatibility of the Abbott IMx Homocysteine Assay with Citrate-Anticoagulated Plasma and Stability of Homocysteine in Citrat Whole Blood. Clin Chem 2001; 47: 1704-1707.
19. Solberg HE. Approved Recommendation (1987) on the Theory of Reference Values. Part 5. Statistical Treatment of Collected Reference Values. Determination of Reference Limits. J Clin Chem Clin Biochem 1987; 25: 645-656
20. Clinical Laboratory Standards/NCCLS. Evaluation of Precision Performance of Clinical Devices- Second Edition; Tentative Guideline. CLSI/NCCLS Document EP5-T2. Villanova, PA: CLSI/NCCLS, March 1992.

EP and US Patents
EP 623174/US 5631127/US 5827645, US 5958717/US6063581

ACL, ACL Futura and ACL TOP are trademarks of Instrumentation Laboratory.

©2006 Instrumentation Laboratory
Issued October 2006

Symbols used / Verwendete Symbole / Símbolos utilizados / Symboles utilisés / Simboli impiegati / Símbolos utilizados / Anvendte symboler / Använda Symboler / Χρησιμοποιηθέντα σύμβολα

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