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1. Aoki N, Moroi M, Yoshida N, Matsuda M. Abnormal Plasminogen. A Hereditary Molecular Abnormality Found in a Patient with Recurrent Thrombosis, J. Clin. Invest. 1978; 61: 1186-1195.
2. Robbins KC. Dysplasminogenemias, Prog. Cardiovasc. Disease. 1992; 4: 295-308.
3. Latallo ZS, Teisseyre E, Lopaciuk S. Evaluation of a Fibrinolytic Profile of Plasma using Chromogenic Substrates. In: Chromogenic Peptide Substrate. Chemistry and Clinical usage. Scully MK, Kakkar VV (eds), Churchill Livingstone, London. 1979; 262-268.
4. Gallimore MJ, Aasen AO, Erichsen NS, Larsbraaten M, Lyngaas K, Amundsen E. Plasminogen Concentrations and Functional Activities and Concentration of Plasmin Inhibitors in Plasma Samples from Normal Subjects and Patients with Septic Shock, Thrombosis Research. 1980; 18: 601-608.
5. Soria J, Soria C, Samama M. A Plasminogen Assay using a Chromogenic Synthetic Substrate. Results from Clinical Work and from Studies of Thrombolytic. In: Progress of Chemical Fibrinolysis and Thrombolysis. Davidson JF, Rowan RM, Samama MM, Desnoyers PC (eds). Raven Press, New York. 1978; 337-346.
6. Richmond JY, McKinney RW eds. Biosafety in Microbiological and Biomedical Laboratories, U.S. Dept. of Health and Human Services, Public Health Service, 4th Edition, 1999.
7. National Committee for Clinical Laboratory Standards. Preparation and Testing of Reagent Water in the Clinical Laboratory, Third Edition, NCCLS Document C3-A3; Vol. 17 No. 18.
8. National Committee for Clinical Laboratory Standards. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays, Fourth Edition, NCCLS Document H21-A4; Vol. 23 No. 35.
9. Zucker S, Cathey MH, West B. Preparation of Quality Control Specimens for Coagulation, Am J. Clin. Pathol. 1970; 53: 924-927.
10. Westgard JO, and Barry PL. Cost-Effective Quality Control: Managing the Quality and Productivity of Analytical Process, AACC Press 1986.
11. IFCC Committee Report, J. Clin. Chem. Clin Biochem. 1983; 21: 749-760.

ACL, ACL FUTURA and ACL TOP are registered trademarks of Instrumentation Laboratory

©1998 Instrumentation Laboratory
Issued July 2006

Symbols used / Verwendete Symbole / Símbolos utilizados / Symboles utilisés / Simboli impiegati / Símbolos utilizados / Anvendte symboler / Använda Symboler / Χρησιμοποιηθέντα σύμβολα

<p>IVD</p> <p><i>In vitro</i> diagnostic medical device <i>In-vitro</i> Diagnostikum De uso diagnóstico <i>in vitro</i> Dispositif médical de diagnostic <i>in vitro</i> Per uso diagnóstico <i>in vitro</i> Dispositivo médico para utilização em diagnóstico <i>in vitro</i> "in vitro" diagnostisk udstyr <i>In vitro</i> diagnostisk medicinsk produkt Προϊόν για διαγνωστική χρήση <i>In vitro</i></p>	<p>LOT</p> <p>Batch code Chargen-Bezeichnung Identificación número de lote Désignation du lot Numero del lotto Número de lote Batch nr. Tillverkningskod Αρ. Παρτίδας</p>	<p></p> <p>Use by Verwendbar bis Caducidad Utilisable jusqu'à Da utilizzare prima del Data límite de utilização Anvendelse Användning Χρήση έως</p>	<p></p> <p>Temperature limitation Festgelegte Temperatur Temperatura de Almacenamiento Températures limites de conservation Limiti di temperatura Límite de temperatura Temperatur begrænsninger Temperatur gräns Περιορισμοί θερμοκρασίας</p>	<p></p> <p>Consult instructions for use Beilage beachten Consultar la metódica Lire le mode d'emploi Vedere istruzioni per l'uso Consultar as instruções de utilização Se vejledning for anvendelse Ta del av instruktionen före användning Συμβουλευτήτε τις οδηγίες χρήσης</p>	<p>CONTROL</p> <p>Control Kontrollen Control Contrôle Controllo Controlo Kontrol Kontroll Υλικό ποιοτικού ελέγχου</p>	<p></p> <p>Biological risks Biologisches Risiko Riesgo biológico Risque biologique Rischio biologico Risco biológico Miljø oplysninger Biologiska risker Βιολογικοί κίνδυνοι</p>	<p></p> <p>Manufacturer Hergestellt von Fabricado por Fabricant Prodotto da Fabricado por Producent Tillverkare Κατασκευαστής</p>	<p>EC REP</p> <p>Authorised representative Bevollmächtigter Representante autorizado Mandataire Rappresentanza autorizzata Representante autorizado Leverandør Auktoriserad representant Εξουσιοδοτημένος αντιπρόσωπος</p>
--	--	--	---	---	--	---	--	---