REF 0501008

CAL 0,5 mL

IVD **CE**

HORIBA ABX
BP 7290 - 34184 Montpellier
cedex 4 - France

Date: 10/03/06

S.A.S au capital de 41.700.000 Euros - RCS Montpellier 328 031 042 - SIRET 328 031 042 000 42 - APE 332 B

Exclusive use:

Minos
Argos & Helios
Micros (8p)
Micros CRP
Micros CRP 200
Pentra 60/60 C+
Pentra 80/XL 80
Pentra 120
Pentra 120 Retic
Pentra DX 120/DF 120
Slide Preparation System

Hematology ABX Devices
(for in vitro diagnostic use)

ABX CRP Cal

EN

1. Functions

ABX CRP Cal is used for calibration of CRP concentration measurement on MICROS CRP instruments.

2. Summary

ABX CRP Cal is a human serum – based Calibrator, used on MICROS – CRP instruments requiring calibration on a periodic basis.

The CRP values are defined by standard IFCC method, by using international reference CRM 470.

3. Principles of measurement

ABX CRP Cal allows calibration of MICROS CRP, previously to measurements taken on whole blood sample (see «Instructions for use» chap 5). The CRP in whole blood sample is measured by immuno-turbidimetry using CRP 100 reagent.

Refer to CRP 100 reagent leaflet for principles of measurement.

4. Warning and precautions

Potentially biohazardous material. For in vitro diagnostic use.

Each human blood donor unit used in the compounding of the serum base for the above referenced products has been tested by and FDA approved method and found to be non reactive for HIV 1/2 (HIV) Antibody, HBc Antibody, HTLV-1 Antibody, HCV Antibody, a Serological test for Syphilis (STS), Hepatitis B Surface Antigen (HBsAg) and HIV-1 Antigen (HIV Ag). The above referenced products were only formed with these donors.

Because no known test method can provide total assurance that products derived from human blood will not transmit infectious diseases, products containing materials from human sources should be handled as if potentially infectious.

Use safe laboratory procedures as outlined in Biosafety in Microbiological and Biomedical Laboratories.

Contain < 0,1% Sodium Azide as preservative. Avoid swallowing and contact with skin or mucous membranes.

5. Instruction for use

See : «User Manual MICROS-CRP», Chapter 4 : Functions / CRP Calibration

6. Storage conditions

Storage temperature: 2 – 8°C

Expiration date: refer to «expiration date» packaging label.^a

Stability after opening: 2 weeks

7. References

Biosafety in Microbiological and Biomedical Laboratories. DC/NIH manual, 1993.

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