REF	0501009: Trol H 0501010: Trol L
CONTROL	0,5 mL
IVD CE	
HORIBA ABX BP 7290 - 34184 Montpellier cedex 4 - France	
Date: 10/03/06	
S.A.S au capital de 41.700.000 Euros - RCS Montpellier 328 031 042 - SIRET 328 031 042 000 42 - APE 332 B	
Exclusive use:	
	Minos Argos & Helios Micros (8p) Micros CRP Micros CRP 200 Pentra 60/60 C+ Pentra 80/XL 80 Pentra 120 Pentra 120 Retic Pentra DX 120/DF 120 Slide Preparation System

A01A00056B

www.horiba-abx.com

EN

1. Functions

ABX CRP Trol L and H are used for accuracy control of CRP concentration measurement on MICROS CRP instruments.

2. Summary

ABX CRP Trol H and ABX CRP Trol L are human serum – based and control fluid, used on MICROS CRP instruments requiring control on a periodic basis.

The CRP values are defined by standard IFCC method, by using international reference CRM 470.

3. Warning and precautions

Potentially biohazardous material.
For in vitro diagnostic use.
Each human blood donor unit used in the compounding of the serum base for the above referenced products has been tested by and FDA approved method and found to be non reactive for HIV 1/2 (HIV) Antibody, HBc Antibody, HTLV-1 Antibody, HCV Antibody, a Serological test for Syphilis (STS), Hepatitis B Surface Antigen (HBsAg) and HIV-1 Antigen (HIV Ag). The above referenced products were only formed with these donors.
Because no known test method can provide total assurance that products derived from human blood will not transmit infectious diseases, products containing materials from human sources should be handled as if potentially infectious.
Use safe laboratory procedures as outlined in Biosafety in Microbiological and Biomedical Laboratories.
Contain < 0,1% Sodium Azide as preservative.
Avoid swallowing and contact with skin or mucous membranes.

4. Instruction for use

See : «User Manual MICROS-CRP», Chapter 4 : Functions / CRP Calibration

5. Storage conditions

Storage temperature : 2 – 8°C

Expiration date: refer to «expiration date» packaging label.^a

Stability after opening :2 weeks

6. References

Biosafety in Microbiological and Biomedical Laboratories. DC/NIH manual, 1993.

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