MS 10310030078
MS 10310030093*

VDRL

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ENGLISH

SUMMARY

Syphilis is a sexually transmitted disease (STD) caused by the spirochete, *Treponema pallidum*.

Other transmission ways are infected blood transfusion, which is rare today due to serological screening, and perinatal (congenital syphilis) that is transmitted from mother to child in utero. Primary syphilis is manifested after an incubation

period of 10-90 days. During the initial incubation period, patients are asymptomatic. The sore, called a chancre, is a firm, painless skin ulceration occurring in 85% of patients. Secondary syphilis is characterized after 4 to 10 weeks after primary infection. The most common symptoms are weight loss, headache, anorexia, myalgia, arthralgia, malaise, fever, lymph nodes and exanthem, occurring in 75 to 100% of patients. On this phase, central nervous system can also be affected. Latent syphilis is after the primary and secondary stages and is characterized by positive serological tests and the absence of clinical findings. This stage can last 1 to 2 years. Tertiary syphilis is manifested in one third of patients without treatment and can present gumma (15%), cardiovascular syphilis (10%) or neuro-syphilis (8-10%).

Serological tests for syphilis are classified as non-treponemal such as **VDRL** (*Veneral Disease Research Laboratory*) and **RPR** (*Rapid Plasma Reagin*) and treponemal tests such as **TPHA** (*Treponema pallidum Hemagglutination*), **FTA-Abs** (*Fluorescent Treponemal Antibody Absorption*) and **ELISA** (*Enzyme-Linked Immunosorbent Assay*).

VDRL from **WAMA Diagnóstica** is a flocculation test, non-treponemal, for syphilis diagnosis by determination of antibodies (reagin) in serum, plasma or CSF. This kit is a modified form of classic **VDRL** since it consists of a stabilized suspension and it is ready for use.

PRINCIPLE OF THE METHOD

When the antigenic suspension of **VDRL** is mixed with serum, plasma or CSF containing antibodies (reagin), the antigen particles will flocculate and the result of the reaction is seen under the microscope. A absence of flocculation indicates a negative result.

KIT PRESENTATION

REF 55250-F (250 determinations)

1. Antigenic suspension (5ml)
2. Instructions for use

REF 55250C-F (250 determinations)

1. Antigenic suspension (5ml)
2. Positive serum control (0.5ml)
3. Negative serum control (1.0ml)

REF 55500-F (500 determinations)

1. Antigenic suspension (10ml)
2. Instructions for use

REF 55500C-F (500 determinations)

1. Antigenic suspension (10ml)
2. Positive serum control (1ml)
3. Negative serum control (1ml)

4. Instructions for use

MATERIAL REQUIRED BUT NOT PROVIDED

- Test tubes
- Serological pipettes
- Tube shelf and rack for tips
- Disposable material container
- Saline 0.9%

STABILITY AND STORAGE

• **ANTIGENIC SUSPENSION (1):** Allow to reach room temperature prior to use. Gently mix prior to use. Stable if stored at 2-8°C up to expiration date. Do not freeze.

• **POSITIVE SERUM CONTROL (2):** Ready for use. Stable if stored at 2-8°C up to expiration date. Contain sodium azide 0.1%. Allow to reach room temperature prior to use.

• **NEGATIVE SERUM CONTROL (3):** Ready for use. Stable if stored at 2-8°C up to expiration date. Sodium azide 0.1%. Allow to reach room temperature prior to use.

The kit presents good performance after begin used for the first time. It is stable up to the expiration date if store at 2-8°C.

SPECIMEN COLLECTION AND STORAGE

Do not use haemolysed, contaminated or lipaemic serum, plasma or LCR for testing. If longer storage is required, store at -20°C from 4 to 6 weeks.

PROCEDURE

a. Qualitative test

Objetivo: It is intended to screen and remove negative serum.

1. Pipette 50µl of the specimens and control sera on the slide ring.
2. Pipette 20µl of shaken antigenic suspension in the same rings. There is no need to mix the components.
3. Rotate the slide test for 4 minutes at 180 r.p.m.
4. Immediately after 4 minutes, inspect the result microscopically.

Obs.: In order to avoid prozone effect we suggest the qualitative test be done with serum, plasma or CSF pure and diluted at 1:8.

READING

Negative: ABSENCE of aggregates. Homogeneous aspect.

Positive weak: PRESENCE of finely dispersed aggregates.

Positive: PRESENCE of medium and large aggregates.

b. Semi-Quantitative test

1. Using a saline solution, prepare serial dilutions of the patient's specimen (1:2, 1:4, 1:16, 1:32 and so on).
2. Pipette 50µl of each dilution on each ring slide.
3. Pipette 20µl of shaken antigenic suspension on each dilution. There is no need to mix these two components.
4. Rotate the test slide for 4 minutes at 180 r.p.m.
5. Immediately after 4 minutes, inspect the result microscopically.

TITER OF SAMPLE

The titer is the last dilution in which the presence of aggregates is still visible.

INTERPRETATION

False positive results can occur with VDRL under such conditions: immunizations, infections, pregnancy, malaria, autoimmune diseases (SLE), malign disease etc. A confirmatory test, specific for treponema, is recommended if VDRL result is positive. WAMA Diagnóstica manufactures and markets FTA-Abs (immunofluorescence) and TPHA (hemagglutination).

TEST PERFORMANCE

VDRL from WAMA Diagnóstica presented satisfactory sensitivity. In 103 known positive samples, the kit showed flocculation from 1+ to 4+. In 146 known negative samples, no flocculation was found.

PRECAUTIONS AND WARNINGS

1. Store the reagents at 2-8°C. Do not freeze.
2. VDRL antigenic suspension contains thimerosal as preservative which may be toxic if ingested.
3. The source material used in the preparation of the controls has been tested and confirmed negative for HCV, HIV and HBsAg. Because no test can offer complete assurance that products derived from human source will not transmit infectious agents it is recommended that these controls be treated as potential biohazards materials.
4. Only for in vitro diagnostic. Do not use the reagents beyond the expiration date.
5. Disposal in accordance with local regulations.
6. Follow the good laboratory practices (GLP) related to storage, handling and material disposal.

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