Imuno-Rápido HBsAg

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ENGLISH

SUMMARY

Hepatitis B is defined as a chronic inflammation of the liver. Most cases of acute viral hepatitis are caused by Hepatitis B Virus (HBV). The complete hepatitis B virion is called the Dane particle. It is DNA virus of the *Hepadnaviridae* family and consists of an outer envelope and core. The outer envelope consists of antigenic proteins and it contains the surface antigen, HBsAg.

Core consists of inner envelope, DNA, and DNA polymerase enzyme. The antigens HBsAg (hepatitis B core antigen) and HBcAg (Hepatitis B e antigen) belong to core.

HBV infection is mainly transmitted by sexual contact, oral sex, perinatal and unsafe injections and transfusions.

Sexual transmission is more frequent with hepatitis A,B,C and Delta virus. Hepatitis B and C are likely to evolve chronic hepatitis and they are also

associated with primary hepatocellular carcinoma.

In Brazil, HBV infection is 8.0% on Amazon, 2.5% in the Southeast and 1.0% in the South.

Acute hepatitis B incubation period is from 45 to 180 days. On symptomatic patients, hepatitis B is likely to evolve to prodromal, icteric and convalescence phases. Hepatitis B is likely to chronically evolve to active or persistent forms, this is demonstrated by laboratorial markers, hepatic function assays and abnormal histology.

Laboratorial diagnostic is performed through serological markers of HBV. Hepatitis B surface antigen (HBsAg) is the first detectable viral antigen to appear during infection, generally it precedes hepatitis clinically evident and also is present in chronic patients.

On the other hand, antibody to the hepatitis B surface antigen (Anti-HBs) can appear on the convalescence phase and it indicates immunity.

The concentration of HBsAg in chronic patients is rarely lower than 20IU/ml.

Imuno-Rápido HBsAg is a test for the qualitative detection of hepatitis B surface antigen HBsAg in serum specimens using a combination of monoclonal and polyclonal antibodies for the selective identification of elevated levels of HBsAg.

PRINCIPLE OF THE METHOD

Hepatitis B surface antigens present on the specimen bind to gammaglobulin conjugate forming an antigen-antibody mixture.

This mixture migrates by capillary action on the test device and binds to anti-HBsAg antibodies in the test region (T). A light pink band appears. The absence of Hepatitis B surface antigen, no colored band will appear in the test region (T). The mixture continues to migrate up to reach control region (C). The conjugate which is not bound to the antigen binds to the reagents and a light pink band appears, showing that the reagents are properly running.

KIT PRESENTATION

REF 617010-R (10 determinations)

1. Test device: 10 units.

2. Instructions for use.

REF 617020-R (20 determinations)

1. Test device: 20 units.

2. Instructions for use.

REF 617040-R (40 determinations)

1. Test device: 40 units.

2. Instructions for use.

MATERIAL REQUIRED BUT NOT PROVIDED

• Paper towel.

• Specimen collection container.

• Automatic pipette.

STABILITY AND STORAGE

• **TEST DEVICE:** It should be kept at 2-30°C. **Do not freeze.**

SPECIMEN COLLECTION AND HANDLING

Only serum specimen should be used with this test. Do not use haemolysed, lipaemic or contaminated serum. Fresh specimen should be used. If longer storage is required, store at 2-8°C for up to 3 days or in freezer at -20°C. Be sure the specimens are brought to room temperature prior to use.

Avoid thawing and freezing repeatedly.


ATTENTION: If the specimen is frozen, it should be brought to room temperature and gently mixed. After this, place it upright for sedimentation.

Do not use specimens of plasma or whole blood, the anticoagulant present in these specimens can interfere on the conjugate action. Diluted sample can cause false negative result.

PROCEDURE

1. Remove the test device (1) from the sealed pouch.

2. Transfer (100µl) of serum to the specimen well of the test device.


3. The result should be read from 10 to 20 minutes. Do not interpret result after 30 minutes.

Obs.: Imuno-Rápido HBsAg detects specimen of serum with hepatitis B virus concentration from 10IU/ml.

READING


NEGATIVE

One light pink band appears in the control region (C).

POSITIVE

Two light pink band appear, one in the test region (T) and the other in the control region (C).

NOTE: The intensity of the pink bands in the control region and in the test region may vary depending on the concentration of HbsAg present in the specimen. This variation does not affect the interpretation of the results.

INVALID

No band appears in the test region (T) and in the control region(C) or no band appears in the control region (C).

EXPECTED VALUES

Immuno-Rápido HbsAg was compared to the electrochemoluminescence method and the results of this comparison is below:

PERFORMANCE CHARACTERISTICS

Sensitivity, Specificity and Clinical Accuracy.

Immuno-Rápido HbsAg identified correctly the specimens from a leading commercial panel and it has been compared with a commercial HbsAg test, electrochemoluminescence method, using clinical specimens. The results showed that the sensitivity of the test is 100% and the specificity is 99%.

Method		Electroquemo-luminescencia		Total Results
Immuno -RÁPIDO	Results	Positive	Negative	
HbsAg	Positive	150	1	151
	Negative	0	99	99
Total results		150	100	250

Sensitivity: 100%

Specificity: 99%

PRECISION

Intra-assay

With-run precision has been determined by testing 10 replicates of four specimens: a negative, a low positive, a medium positive and a high positive. All specimens were correctly identified in more than 99% of the time.

Inter-assay

Between-run precision has been determined by 10 independent assays on the same four specimens: a negative, a low positive, a medium positive and a high positive. Three different lots of **Immuno-Rápido HbsAg** have been tested over a 3-month period using a negative, a low positive, a medium positive and a high positive specimens. All specimens were correctly identified more than 99% of the time.

LIMITATION

Immuno-Rápido HbsAg is a screening test for the detection of hepatitis B surface antigen, therefore, a repeated positive result is an evidence of this virus on specimen. Positive results should be confirmed with other test.

False positive and false negative results can occur with this test. The proportion will depend of the prevalence of the disease on the tested population.

PRECAUTIONS AND WARNINGS

- For professional "in vitro" diagnostic use only.
- Read carefully this instructions for use insert prior to testing.
- The expiration date is the last day of the month printed on test device pouch and on the kit box.
- Keep the kit away from heat and sunlight.
- Do not freeze the test device. Irreversible deterioration may be yielded.
- Sodium azide 0.1% may react with lead and cooper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide might be toxic if ingested.
- Allow the kit components to equilibrate to room temperature prior to testing.
- Do not use the kit beyond the expiration date.
- Do not interchange kit components with those from sources other than the same catalog number from WAMA DIAGNOSTICA.
- Disposal in accordance with local regulations.
- Follow the good laboratory practices (GLP) in storing, dispensing and disposing of these materials.

WARRANTY

WAMA Diagnóstica replaces this kit since it is not beyond expiration date. The

returned kit must be evaluated by Wama's technical support. The warranty will be invalid and kit will not be replaced if technical support finds evidence that running, handling and storage were not properly followed.

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