Imuno-Rápido HCG

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ENGLISH

SUMMARY

The determination of human Chorionic Gonadotropin (hCG) is universally recognised as diagnostic test for pregnancy.

Human Chorionic Gonadotropin is a glycoprotein hormone, molecular weight around 37,000 daltons, secreted by trophoblastic cells of the placenta during pregnancy.

Chorionic gonadotropin is made of 2 different chains known as alpha and beta. Alpha is identical to LH molecule whereas beta differentiates from alpha due to 30 aminoacids, present in the carboxy terminal but absent in the LH molecule.

hCG is secreted from 6 to 8 days after conception and the levels continue to rise very rapidly peaking in the 50,000-200,000mIU/ml range about 6th to 8th week and then stabilising to an average level of 5,000-20,000mIU/ml up to the end of pregnancy.

Imuno-RÁPIDO HCG is one step chromatographic immunoassay test to selectively detect HCG on specimen with sensitivity of 25mIU/ml.

PRINCIPLE OF THE METHOD

HCG present on specimen binds to monoclonal antibody-anti- β -hCG-colloid gold conjugate forming an antigen-antibody complex.

The mixture then migrates on the membrane of the test strip/device and binds to antibodies anti-hCG in the positive reaction region (T) producing a light-pink colored band if the concentration is greater than 25mIU/ml. The absence of this colored band is indicative of a negative result. The mixture continues to migrate to control region (C). The conjugate which is not bound to antigen joins to the reagents producing a light-pink colored band in the control region (C) serving as a procedural control and indicating that the test has been performed properly.

KIT PRESENTATION

REF 620040T-R (40 determinations)

1. Test-strip: 40 units

2. Instructions for use.

REF 620080T-R (80 determinations)

1. Test-device: 80 units

2. Instructions for use.

REF 620010-R (10 determinations)

1. Test-strip: 10 units

2. Instructions for use.

REF 620020-R (20 determinations)

1. Test-strip: 20 units

2. Instructions for use.

REF 620040-R (40 determinations)

1. Test-device: 40 units

2. Instructions for use.

MATERIAL REQUIRED BUT NOT SUPPLIED

- Micropipettes.

- Paper towel.

- Test tube and rack.

- Disposable material container.

STABILITY AND STORAGE

Test-strip/Test-device: Keep at 15° to 30°C. Do not freeze.

SPECIMEN COLLECTION AND STORAGE

Tests can be performed in serum or urine.

Serum

Serum free of lipaemy, haemolysis and contamination are recommended.

Keep the serum at 2° to 8°C for 48 hours. For longer storage, keep it in the freezer at -20°C for up to 1 year.

Frozen specimens should be mixed and thawed at room temperature prior to use. Avoid freezing and thawing repeatedly and foam formation.

Urine

The specimens of urine must be collected in a clean container, without preservative. Specimens collected at any time may be used. However, the first morning urine generally contains the highest concentration of hormone.

The specimens not tested after the collection should be stored refrigerated at 2° to 8°C for 24 hours. If longer storage is required, the specimens should be stored in freezer. Frozen specimens should be mixed and thawed at room temperature prior to use. Avoid freezing and thawing repeatedly and foam formation.

PROCEDURE


Test-strip

1. Dip the strip into the specimen

for 10 seconds.

Avoid surpassing the limit indicated

by arrows.

2. Wait for 5 minutes. Read it.


Obs.: Reading after 10 minutes is not recommended.

Test-device

1. Pipette 100 µl of serum or urine (avoid


air bubbles) in the specimen

well (►) of test-device.


2. Wait for 5 minutes. Read it.
Obs.: Reading after 10 minutes is not recommended.


READING
NEGATIVE: Only one colored band appears in the control region (C).

POSITIVE: It will appear two colored bands, one in the test region (T) and other in the control region (C).

INVALID: If no band appears in the test (T) and control (C) regions, or if no band appears in the control region (C).

TEST PERFORMANCE

Specificity:
 Cross-reaction on the following concentrations of hormones do not interfere Imuno-Rápido HCG from Wama.

hTSH	1.000mIU/ml	WHO 6838
hLH	500mIU/ml	WHO 2° IRP HMG
hFSH	1.000mIU/ml	WHO 2° IRP HMG

Postmenopausal urine: A study performed using a total of 35 randomly selected specimens of urine from postmenopausal women showed negative results in all specimens. This study was done since the urine from women in menopause interferes on pregnancy tests due to cross reaction with gonadotropic hormones.

Interference testing: The following substances were added to specimens of urine with hCG levels of 0 and 50mIU/ml. No changing on the expected results was indicated with substances and concentrations below:

Acetaminophen 20mg/dl	Caffeine: 20mg/dl	Ampicillin: 20mg/dl
AAS: 20mg/dl	Gentiesic Ac.: 20mg/dl	Tetracilin: 20mg/dl
Ascorbic Acid: 20mg/dl	Glucose: 2g/dl	
Atropine: 20mg/dl	Hemoglobin: 1mg/dl	

Precision:

Intra-assay: Within-run precision was determined by 10 replicates of three specimens with hCG ranging from 0.20 to 100mIU/ml. The negative and positive values were correctly identified in 100% of the time.

Inter-assay: Between-run precision was determined by 10 independent assays with three different lots of Test-strip/Test-device for 10 months using the same three samples ranging from 0.20 to 100mIU/ml of hCG. The negative and positive values were correctly identified in 100% of the time.

Sensitivity:

Imuno-Rápido HCG detects levels of hCG of 25mIU/ml in accordance with WHO International Standard.

Occasionally, specimens of urine or serum with levels lower than 25mIU/ml can show positive results.

Specimens with high levels of hCG (200,000 mIU/ml) show positive results when they are tested.

Urine from non-pregnant women and healthy men normally show no detectable levels of hCG when tested with **Imuno-Rápido HCG**.

Although the levels of hCG vary greatly with gestational age and between individuals, Imuno-Rápido HCG from WAMA is capable of detecting pregnancy as early as 1 day after the first missed menses.

USE LIMITATION

This test does not distinguish normal from ectopic pregnancy. A number of conditions other than pregnancy, including trophoblastic disease and certain non-trophoblastic neoplasms which includes choriocarcinoma and hydatid mole, cause elevated levels of hCG. Therefore, the presence of hCG in urine specimen should not be used to diagnose pregnancy unless these conditions have been ruled out. Negative results are expected at the low gestational age, which may not contain representative levels of hCG. On this situation, a first morning urine specimen should be collected 48 hours later and tested. High levels of hCG can be detected for some weeks after natural childbirth, caesarean, and therapeutic or spontaneous abortion.

Some specimens of serum with high concentration of rheumatoid factor can show positive results non-specific when tested.

Specimens of urine which are not collected from the first morning urine may not contain representative levels of hCG for testing. If the pregnancy is still suspected

even with a negative result, a first morning urine specimen should be collected 48 hours later and tested. As with all diagnostic procedure, the result of this test should be interpreted with other clinical information available to the physician.

PRECAUTIONS

- For in vitro diagnostic use only for use.
- Read carefully the instructions prior to testing
- Protect the kit from high temperatures and direct sunlight
- Do not freeze the Test-strip/Test-device.
- Store the kit at 15° to 30°C. Do not store it in the refrigerator.
- Do not use the kit after expiration date.
- Disposal material in accordance with local regulations.
- Follow good manufacturing practice (GMP) on storing, handling and material disposable.

WARRANTY

WAMA Diagnóstica replaces this kit since it is not beyond expiration date. The returned kit must be evaluated by Wama's technical support. The warranty will be invalid and kit will not be replaced if technical support finds evidence that running, handling and storage were not properly followed.

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