Imuno-Rápido HIV

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ENGLISH

SUMMARY

Since the acquired immunodeficiency syndrome, AIDS, first appear in 1981, the human immunodeficiency virus (HIV) infection became a global epidemic. Although intravenous drug users, unprotected sexual intercourse, infected mother to her baby and through blood transmission, considered as a risk group, unprotected sex (homosexual and are heterosexual) is responsible for the majority of infection cases among adults around the world. HIV is a viral RNA genome from Lentiviridae family and belongs to non-cytopathic and non-oncogenic retroviruses group. In order to multiply itself, the virus need an enzyme called reverse transcriptase, which is responsible for viral RNA transcription to a DNA copy to the genome host.

HIV-1 virion is cilinder and surrounded by a lipidic envelope. Its RNA genome has 10,000 base pairs and characterize itself by the presence of two LTR (Long Terminal Repeats) and nine regions, three codifier of proteins (Gag, pol and env) and six regulatory (tat, rev, nef, vpr, vif, vif). Gag, pol and env regions codify the proteins and glycoproteins p-24, p17, gp-120, gp41 and the reverse transcriptase enzymes, protease and integrase. There are two types of HIV, known as HIV-1 and HIV-2. The genetic homology seen on proviral genomes is roughly 40-45%, which determine the codification and synthesis of several similar gene products. Such detail justifies a serum cross-reactivity among both, seen on diagnostic serologic tests. Variant of genome (subtypes) of HIV-1&2 have been described. HIV-1 subtypes are classified in two groups: M (major) and O (outlier) with genetic variability up to 30% on env segment. In the M group, there are nine subtypes (A,B,C,D,E,F,G,H and I) and in the O group, there is only one. In relation to HIV-2, there are five subtypes (A,B,C,D and E).

Imuno-Rápido HIV-1 and 2 uses as capture reagent, immobilized on test device, HIV-1 recombinant proteins, which correspond to immuno-dominant regions of gp-41 and an epitope which correspond to gp-36 protein from HIV-2.

PRINCIPLE OF THE METHOD

Anti-HIV-1 and/or anti-HIV-2 antibodies present on the sample bind to conjugate antibody forming a complex. This complex migrates through the membrane and bind to the antigen (gp-41/gp36) in the test region (T), one light pink colored band appears. On the absence of anti-HIV-1 an/or anti-HIV-2 antibodies, no colored band appears. A control reagent immobilized on membrane of test device will determine the appearance of a second colored band showing the reagents are running properly (control region "C").

KIT PRESENTATION

REF 622010-R (10 determinations)

1. Test device: 10 units

2. Diluent solution: 1 x 2ml

3. Instructions for use

REF 622020-R (20 determinations)

1. Test device: 20 units

2. Diluent solution: 2 x 2ml

3. Instructions for use

REF 622040-R (40 determinations)

1. Test device: 40 units

2. Diluent solution: 4 x 2ml

3. Instructions for use

MATERIAL REQUIRED BUT NOT PROVIDED

Paper towel

Disposable material container

Automatic pipette

STABILITY AND STORAGE

TEST DEVICE (1): Keep at 2-30°C. **Do not freeze.**

DILUENT SOLUTION (2): Ready for use. Sodium azide 0,1% as preservative. Stable if stored at 2-30°C up to expiration date. **Do not freeze.** Allow the components to reach room temperature prior to use.

SPECIMEN COLLECTION AND STORAGE

Specimens of serum, plasma or whole blood free of haemolysis, lipaemy and contamination are recommended. Flesh specimens may be stored at 2-8°C for up to 48 hours if not used immediately. For prolonged storage, the specimens of serum and plasma should be frozen at -20°C. Do not freeze whole blood specimens.


ATTENTION: When the specimen is frozen, it should reach room temperature and homogenized completely. Keep the specimen vial

in vertical position to allow any particle to be sedimentated. Do not shake the specimen.

Diluted sample can cause false negative result.

PROCEDURE


1. Remove the test device (1) from the envelope.
2. Pipette 10µl of sample (with no air bubbles) on sample well (▶) of test device.


3. Dispense 3 drops (100µl) of diluent solution (2) in the sample well of test device.
4. Read from 10 to 15 minutes. Do not read after 20 minutes.

READING



NEGATIVE: Appearance of one light pink colored band in the control region (C).

POSITIVE: Appearance of two colored bands, one in the test region (T) and other in the control region (C).

INVALID: No colored bands appear in the test region (T) and in the control region (C), or if no colored band appears in the control region (C).

OBS.: Disregard the results after the prescribed time for the readings.

TEST PERFORMANCE

The **Imuno-Rápido HIV** presented 100% of sensitivity. 137 known positive specimens were used in the quality tests. In 1131 known negative specimens, the Imuno Rápido HIV presented only one sample with false positive result which provided specificity of 99.9%

USE LIMITATION

The **Imuno-Rápido HIV** is a screening test aimed to detect anti HIV-1 and/or HIV-2 antibodies, therefore, a repeatedly positive result is a strong evidence of antibodies on the specimen.

Positive results should be confirmed by Western Blot Method or other confirmatory methods.

False positive and false negative results may occur, the proportion depends on the prevalence of the disease on the tested population. As any other diagnostic procedure, test result must be interpreted with further clinical data.

PRECAUTIONS AND WARNINGS

1. For "in vitro" diagnostic use only.
2. Read carefully the instructions for use prior to use.
3. Expiration date is the last day of the month printed on labels.
4. The test kit should be kept away from direct sunlight
5. Do not freeze test device.
6. The kit contains 0.1% sodium azide as a preservative. Sodium azide may react with lead and copper plumbing to form highly explosive salts. On disposal, flush with large quantities of water. Sodium azide may be toxic if ingested.
7. Do not use kit beyond the expiration date.
8. Do not interchange kit components from different kits.
9. Disposal in accordance with local regulations.
10. Follow good laboratory practices (GLP) related to storage, handling and material disposal.

WARRANTY

WAMA Diagnóstica replaces this kit since it is not beyond expiration date. The returned kit must be evaluated by Wama's technical support. The warranty will be invalid and kit will not be replaced if technical support finds evidence that running, handling and storage were not properly followed.

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